



B.IV4



Strip-type biosensor




For various bioanalytical applications



Benefits and characteristics




- Enzyme-based amperometric measurement
- Excellent long-term stability
- Reference, counter, and blank electrodes on-chip
- Outstanding reliability
- Fast response time
- Gamma and beta sterilization compatible
- Small size
- Suitable for dip-in applications
- For industrial applications



Left: B.IV4 sensor with JST SURS; **Center:** B.IV4 sensor FPC area; **Right:** B.IV4 sensor sensing area


Illustration

B.IV4 sensor, flex print cable termination



B.IV4 sensor with FPC, all dimensions in mm

B.IV4 sensor, JST SURS 4-pin connector



B.IV4 sensor with SURS 4-pin, 0.8 mm pitch connector, all dimensions in mm



Technical data



Dimensions:	See illustration above. Dimensions in mm	
Measurement principle:	Enzymatic amperometric (oxidase enzymes and H ₂ O ₂ oxidation)	
Working electrode:	Platinum covered with enzyme membrane	
Blank electrode:	For background compensation	
Reference electrode:	Silver/silver chloride	
Counter electrode:	Platinum	
Analytes and measurement range:	Glucose	0.1 mM to 50 mM
	Lactate	0.05 mM to 25 mM
	Glutamine	Contact iST Jobst
	Glutamate	Contact iST Jobst
	Notes:	- measured in acetate buffer at 37 °C - other measurement ranges and analytes on request

Sensitivity at 37 °C:	Glucose	Typically: 0.4 nA/mM
Sterilization	Irradiation (beta, gamma)	
	- recommended dose: < 25 kGy	
	- sensitivity increases according to applied dose	
	- lifetime decreases according to applied dose	
	Initial bioburden < 1 CFU per sensor	
	Do not use organic solvents, for more information, contact iST Jobst	

Time (t_{90%}) to first measurement: ~ 30 min after storage at 37 °C

Response time (t_{90%}): < 90 s in acetate buffer at 37 °C

Temperature influence:

Glucose:	~3.8 %/K
Lactate:	~3.2 %/K

Storage conditions: +4 °C to +35 °C, desiccated

Shelf life at recommended storage conditions:

> 6 months (from delivery)
> 2 years from fabrication

Operational lifetime:

Glucose:	> 120 days at 20 mM > 50 days at 50 mM
Lactate:	> 7 days at 8 mM
Notes:	- measured in acetate buffer at 37 °C - lifetime may vary in other buffer systems

Operating temperature: 15-42 °C

pH range: 6-8

Drift at 37 °C: < 5 %/day

Suitable pH buffer systems:

Bicarbonate, acetate, imidazole, for more information, contact iST Jobst

- to be used in buffered media only
- buffer must contain chloride, [Cl⁻] ca. 110 mM
- not suitable for direct use in tap water or DI-water


Electrical connection:

FPC: compatible with 0.5 mm pitch FFC/FPC ZIF connectors
 JST: JST SURS header 4-pin 0.8 mm female


General note: Performance data in this document was determined in acetate buffer at 37 °C, pH 7 and normal atmospheric conditions. All parameters may vary in other media.



Pin assignment



B.IV4 FPC sensor contacts



B.IV4 sensor JST SURS 4-pin 0.8 mm pitch connector side

1	2	3	4
Reference	Counter	Blank	Glucose or Lactate

Order information

Product description	Product name	Product number	Order code iST
Biosensor glucose	B.IV4.G.1X.FPC	1.00102.000	105127
Biosensor glucose, JST connector	B.IV4.G.1X.JST	1.00102.002	105192
Biosensor lactate, JST connector	B.IV4.L.1X.JST	1.00102.003	105191
Other analytes		On request	On request

Disclaimer

Evaluation product for professionals to be used solely for research and development purposes! Not for medical and diagnostic use. Not to be used on humans. For more information, contact iST Jobst.



Jobst Technologies GmbH • Georges-Köhler-Allee 302 • 79110 Freiburg • Germany
+49 761 55 77 520 • info.jobst@ist-ag.com • www.jobst-technologies.com

ALL MECHANICAL DIMENSIONS ARE VALID AT 25 °C AMBIENT TEMPERATURE, IF NOT DIFFERENTLY INDICATED • ALL DATA EXCEPT THE MECHANICAL DIMENSIONS ONLY HAVE INFORMATION PURPOSES AND ARE NOT TO BE UNDERSTOOD AS ASSURED CHARACTERISTICS • TECHNICAL CHANGES OR PRODUCT SPECIFICATIONS WITHOUT PREVIOUS ANNOUNCEMENT RESERVED • THE INFORMATION ON THIS DATA SHEET WAS EXAMINED CAREFULLY AND WILL BE ACCEPTED AS CORRECT; NO LIABILITY IN CASE OF MISTAKES • LOAD WITH EXTREME VALUES DURING A LONGER PERIOD CAN AFFECT THE RELIABILITY • THE MATERIAL CONTAINED HEREIN MAY NOT BE REPRODUCED, ADAPTED, MERGED, TRANSLATED, STORED, OR USED WITHOUT THE PRIOR WRITTEN CONSENT OF THE COPYRIGHT OWNER • ALL RIGHTS RESERVED.